


Product Reference Guide

Glass

Standard Glass Finished Single Vision

70mm Finished Single Vision (White)

Sphere Power Up to -2.00 cyl	Base Curve	Curve Radius	Center Thickness	Sphere Power Up to -2.00 cyl	Base Curve	Curve Radius	Center Thickness
Finished Single Vision, 70mm, Minus Powers				Finished Single Vision, 70mm, Plus Powers			
Plano	6.00	106.000	2.2	0.25	6.00	106.000	2.2
-0.25	5.00	106.000	2.2	0.50	6.00	106.000	2.4
-0.50	5.00	106.000	2.2	0.75	6.00	106.000	2.6
-0.75	5.00	106.000	2.2	1.00	6.00	106.000	2.8
-1.00	5.00	106.000	2.2	1.25	6.00	88.333	3.0
-1.25	5.00	106.000	2.2	1.50	6.00	88.333	3.2
-1.50	5.00	106.000	2.2	1.75	6.00	88.333	3.4
-1.75	4.00	132.500	2.2	2.00	6.00	88.333	3.6
-2.00	4.00	132.500	2.2	2.25	7.00	75.714	3.9
-2.25	4.00	132.500	2.2	2.50	7.00	75.714	4.2
-2.50	4.00	132.500	2.2				
-2.75	4.00	132.500	2.2				
-3.00	4.00	132.500	2.2				
-3.25	3.00	132.500	2.2				
-3.50	3.00	132.500	2.2				
-3.75	3.00	176.667	2.2				
-4.00	3.00	176.667	2.2				
-4.25	3.00	176.667	2.2				
-4.50	3.00	176.667	2.2				
-4.75	3.00	176.667	2.2				
-5.00	2.00	176.667	2.2				



Note: -4.25 to -5.00 available only in 0.00 cylinder.