

Product Reference Guide

Glass

High Index Semi-Finished Single Vision

66mm 1.70 Index (White)

| Base Curve | True Curve (1.530) | Front Radius | 50mm SAG* | Back Curve | Back Radius | CT (mm) Nominal | ET (mm) Nominal |
|------------|--------------------|--------------|-----------|------------|-------------|-----------------|-----------------|
| Plano | 0.00 | - | - | -10.00 | -53.00 | 4.5 | 15.6 |
| 0.50 | 0.50 | 1060.000 | 0.295 | -10.00 | -53.00 | 5.5 | 15.6 |
| 2.00 | 2.00 | 265.000 | 1.182 | -9.00 | -58.89 | 5.3 | 12.4 |
| 4.25 | 4.15 | 127.711 | 2.471 | -8.00 | -66.25 | 5.2 | 8.7 |

*Utilized for power calculation only, not for thickness calculation

Product Reference Guide

Glass

High Index Semi-Finished Single Vision

71mm 1.80 Index (White)

| Base Curve | True Curve (1.530) | Front Radius | 50mm SAG* | Back Curve | Back Radius | CT (mm) Nominal | ET (mm) Nominal |
|------------|--------------------|--------------|-----------|------------|-------------|-----------------|-----------------|
| 0.50 | 0.50 | 1060.000 | 0.295 | -8.00 | -66.25 | 7.8 | 11.0 |
| 2.00 | 2.00 | 265.000 | 1.182 | -8.00 | -66.25 | 6.0 | 13.9 |
| 4.25 | 4.15 | 127.711 | 2.471 | -8.00 | -66.25 | 6.0 | 11.2 |

*Utilized for power calculation only, not for thickness calculation