

# How To Install (DVI RxWizard)

## Prerequisites

---

1. An Active Internet Connection
2. Administrator Rights on the computer you are installing the software


## Steps

---


1. Open Internet Explorer
2. Enter <http://thedvi.com/plugin/> into the address bar and hit "Enter"
3. Click the Link to download the latest release of RxWizard




4. When Prompted, Click "Run"




5. Wait for the program to download
6. If Prompted, Click "Yes"




7. Click "Next"




8. Enter the password and click "Next" (password: slaboff)




9. Click "Next"




10. Click "Yes"




11. Click “Next”



12. Click “Next”



13. Click "Finish"



14. Double-Click on "DVI RxWizard" Desktop Icon



15. Verify Settings are correct and Click "OK"

DVI Remote Entry Options

Default Edge Type

Default Frame

Job Ship Options

Days to keep archived jobs

Autoprint jobs

Autoprint Frame Packing Slip

Transmission Method

Auto-generate Rxnumber  Next

No Frame Database Warning

GeoTrace Digitizer

Prompt for AR Coating

Network does not support file locking

Disable Auto-tabbing

Operator Entry  Operator Setup

Use Curl for Internet Transfers

Account-Selective Transmit

Enhanced Order Lookup   
(Requires data)

Debug (DVI Internal Use)

OK Cancel